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ENGLISH

SURGICAL INSTRUMENTS AND ACCESSORIES

For dental use

INSTRUCTIONS FOR USE

1. GENERAL INFORMATION

All surgical instruments and accessories are made of stainless steel and are invasive, reusable and supplied non-sterile. In order to use the devices correctly, it is necessary to carefully read and apply the Information for the User provided by Prodent Italia. This information includes these Instructions for Use and the relevant clinical/surgical protocol for each device, as indicated in the INDICATIONS paragraph below. The Information for the User can be found in the company website or by contacting Prodent Italia S.r.l. directly.

2. INDICATIONS

Surgical instruments and accessories are intended for use when preparing the surgical site and inserting Prodent Italia endosseous devices in the implant site. The devices can either be specific to a certain implant range, to a certain implant-restoration connection or common to several ranges or connections. Each device must be used in accordance with the surgical procedures described in detail in the clinical/surgical protocol for the devices to be implanted. The specific indications for use are given below and the user should therefore consult the protocol for the implants to be placed: CL008 for PRIME and TWINNER implants, CL010 for CONOMET implants, CL011 for PROSHAPE implants, CL007 for PROVIDE implants and CL012 for PROFAST implants.

CE 0546 devices:

Driver: available as a single variant for each restoration range and for each connection; to be used for picking up the implants, positioning them in the receptor site and starting treatment. Insert the driver connected to the digital wrench in the implant connection. Check that the driver has completely engaged the implant connection; if necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring marked on the driver as reference). Move the implant to the receptor site, position it so that it is firm and stable and disconnect the digital wrench. Connect the driver to the Ratchet and tighten the implant until it is correctly positioned. Follow the torque indications provided in the protocols and Instructions for Use for the implants to be used. The drivers are fitted with O-rings and a dedicated square connector for connection with the digital wrench and the Ratchet. To identify the size variants, refer to the specific protocol for the implant being inserted.

Contra-angle Driver: available in different variants for different restoration ranges; to be used connected to the contra-angle to pick up the implants, position them in the receptor site and insert them without EVER exceeding 25 rpm, and 45 Ncm; Insert the driver connected to the contra-angle in the implant connection. Check that the driver has completely engaged the implant connection; if necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring

marked on the driver as reference). Take the implant to the receptor site and tighten the implant until it is correctly positioned. To identify the size variants, refer to the specific protocol for the implant being inserted.

Extension: available as a single variant; to be used when it is deemed necessary to increase the length of the connection between the digital wrench or the Reverse Ratchet or connection for the contra-angle and dedicated devices. It is fitted with an O-ring for connection with the above devices.

Bur extension: available as a single variant (Ref. Clinical Protocols); to be used when it is deemed necessary to extend the connection between the bur and the contra-angle. Do not use the device beyond 45 Ncm.

Contra-angle connection: available as a single variant; to be used for mechanical use of the screw traps and for mechanical implant insertion. NEVER exceed 25 rpm and 45 Ncm during use.

Devices:

Digital wrench: available as a single variant; it MUST be connected to Drivers, Hex Wrenches and screw traps so that they can be used manually in a correct manner; All the above devices are fitted with O-rings and a dedicated square connector for connection with the wrench.

Hex screwdriver: available in three different lengths to satisfy different anatomical needs; to be used to tighten or loosen the screws as indicated in the relevant protocol; it has a loop hole so that it can be tied to a safety thread. **N.B.:** The devices have two laser-etched notches on the shaft between the tip and the handle, in order to distinguish them from the screwdrivers of previous version, which do not have this kind of mark. **N.B.:** **DO NOT** use screwdrivers without notches to tighten devices with TS connections (for CONOMET and PROSHAPE implants)

Screwdriver Premilid: longer variant than hex screwdrivers; suitable for tightening or loosening intact and customized Premilid abutments.

Extraoral wrench: available as a single variant; to be used with the Hex Screwdrivers in order to keep them in the correct position during tightening; the pin perpendicular to the body of the instrument must be inserted in the hole on the head of the screwdriver.

3. CONTRAINDICATIONS

All surgical instruments and accessories must be used exclusively according to the indications provided and to prepare the site and insert the devices they are to be used for. The contraindications concerning Prodent Italia implantable devices are listed in the corresponding Information for the User.

4. SIDE EFFECTS AND COMPLICATIONS

After oral surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. The long-term risks that can be associated to implant-restoration or guided bone regeneration include bone reabsorption, failure of integration, oedema, chronic pain and dehiscence. The success of implant-restoration or guided bone regeneration treatment is closely related to the correctness of all the phases involved, including the preparation of the receptor site. Correct and careful use of surgical instruments and accessories is a key factor for reducing the risks of implant-restoration or guided bone regeneration treatment failure, in addition to reducing the severity of temporary side effects. Practitioners are asked to report any cases of failure to Prodent Italia S.r.l

5. PACKAGING AND STERILITY

All devices are provided NON-STERILE and individually packed; depending on the nature of each one, they may be supplied in a plastic container inside sealed secondary packaging, or directly inside heat-sealed transparent plastic bags or directly inside a sealed box. The label identifying each device is placed, according to the cases, on the secondary packaging, box or plastic bag. No packaging is suitable for sterilisation.; before sterilising the devices, remove them from the packaging and clean them, as indicated in the following paragraph.

6. CLEANING AND STERILISATION

All devices are reusable. Immediately before use, they MUST be **cleaned and sterilised, using a valid method and under the practitioner’s responsibility** (Prodent Italia recommends compliance with standard UNI EN ISO 17665 for the development of validation and routine control of the steam sterilisation process, and with standard UNI EN 13060 for determining the testing methods for steam sterilising equipment and for defining the sterilisation cycles to be used). Immediately after use, the devices MUST be cleaned to remove all types of residue, using plastic brushes with stiff, non-metal bristles. We recommend washing in an ultrasonic bath using appropriate detergents, avoiding using hydrogen peroxide, oxidising acids (such as oxalic, sulphuric and nitric acids) and high concentrations of chlorine, which can damage the devices. Prodent Italia S.r.l. provides further information on the washing and sterilising processes in the *“Important recommendations for the cleaning and sterilisation of devices”* section of the relevant protocol. Prodent Italia S.r.l. has tested the efficacy, on devices of its own manufacture, of the sterilisation process in a saturated steam autoclave at 134°C for 5 minutes.

⚠️ WARNINGS

- Before using surgical instruments and accessories read the updated Information for the User carefully.
- The name on the device label alone may not be sufficient for identifying its intended use: see Information for the User.
- All devices MUST ONLY be used for dental surgical treatment, by qualified dental professionals who are fully knowledgeable with the theoretical and practical aspects of surgical procedures, and under their own responsibility.
- All devices must only be used according to the indications provided and with the devices they are intended for. When available as different sized variants, identify using markings and colour-coding.
- NEVER use devices whit marking when the markings are deteriorated.
- ALWAYS check that the packaging is undamaged before using the device the first time. NEVER use devices with visibly damaged packaging.
- ALWAYS check the device is undamaged before each use. NEVER use devices that are visibly damaged. If in doubt as to the integrity of the devices, request a check by Prodent Italia.
- Take the greatest care during the washing and sterilisation phases, which MUST be performed using a validated method by and under the responsibility of the practitioner.
- NEVER use devices fitted with an O-ring, if it is missing or worn.
- Before use, it is the practitioner’s responsibility to make sure that each device is fit for use and suitable for the intended use and for the envisaged purposes, also according to the Information for the User.
- When connecting two or more devices, ALWAYS check that they are correctly and fully connected before and during use as specified in the indications for use.
- We recommend securing the instruments, which have a loop hole, with a safety thread.
- All devices should be handled with special care, to avoid accidental damage that could cause them not to function properly.
- Handle all devices in aseptic conditions and wearing sterile gloves.
- NEVER force tightening.
- The practitioner is responsible for planning the treatment, choosing the surgical technique, the type and size of the device to be implanted, the instrument set and its use, the rehabilitation technique and for the planning and performing of routine checks. Mistakes and inaccuracy during these phases may lead to a failure of the implant-restoration treatment.
- The device MUST be handled as hospital waste in the event of disposal. The device should therefore only be disposed of after regular cleaning, in a controlled manner and according to current applicable legislation, whenever necessary.
- Prodent Italia S.r.l. ONLY accepts returned devices in their undamaged original packaging and within the contractual terms provided.
- Prodent Italia S.r.l. declines all responsibility for failure to comply with the updates to the Information for the User.
- Prodent Italia S.r.l. furthermore declines all express or implicit liability in relation to direct, indirect or any other type of damage which may be connected to or arise from any evaluation or professional practice error, made during the use of the various products manufactured by Prodent Italia. The use of dental implants, all surgical treatments and prosthetic restorations must always be made under the practitioner’s control, who therefore assumes full liability thereof. The user of Prodent Italia’s items is solely responsible for evaluating whether a product is suitable or not for each individual patient and clinical case.

KEY for the symbols present on the label:

LOT	Batch code.
REF	Catalogue number.
⚠️	Caution! The device is associated with warnings not indicated on the label and that are described in the Instructions for Use.
ℹ️	Consult the Instructions for Use inside the pack
🚫	DO NOT use if the package is damaged
🏭	Manufacturer pursuant to 93/42/EEC and subsequent amendments

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