Via Pitagora, 9 – 20016 Pero (MI) Italy – www.prodentitalia.eu

ITALIANO

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ENGLISH

PROSTHETIC COMPONENTS AND FASTENING SCREWS

For dental use

INSTRUCTIONS FOR USE

1. GENERAL INFORMATION

CONIK-FIT prosthetic components are SINGLE USE, supplied non-sterile and available in various types and sizes, to make prostheses that meet the different anatomical and functional requirements. CONIK-FIT prosthetic components are divided into reformed prosthetic components, which are invasive, and calcinable components, which are laboratory devices. The devices to be connected to the implant using fastening screws are supplied together with the dedicated screw. The table lists the types of prosthetic components, the materials they are made of and the dedicated accessories they are supplied with.

Restoration components	Dedicated accessories (in packaging)	Material
PREFORMED  0546		
Straight abutments	Fastening screw	Grade 5 titanium
17° Abutments		
Milling abutments		
Premilled abutments		
Multi abutment		
Aesthetic bases		
Temporary aesthetic abutments	Fastening screw	Grade 5 titanium base PEEK bush
Temporary abutment peek sleeve	/	PEEK
CASTABLES (laboratory devices) 		
Castable abutments	Fastening screw	
Cementable castable abutments	/	Polycarbonato
Multi abutment castable sleeves	/	

All fastening screws for CONIK-FIT prosthetic components are made of grade 5 titanium, they are SINGLE USE and supplied non-sterile. All the devices are part of the CONIK-FIT implant-restoration system. To use them correctly, carefully read and apply the Information for the User provided by Prodent Italia. This information consists in the Instructions for Use and the CL007 Clinical Protocol, all of which can be found in the company website or by contacting Prodent Italia S.r.l. directly. The protocol indicates the types and sizes of the CONIK-FIT prosthetic components and the corresponding fastening screws, needed to meet the different anatomical and restorative requirements of each of the CONIK-FIT prosthetic ranges. CONIK-FIT prosthetic components are available only for one prosthetic range, Q3,5 common to all PROVIDE implant diameters, and DO NOT have markings or colouring.

2. INDICAZIONI

CONIK-FIT prosthetic components are interconnection devices between the endosseous implants of the CONIK-FIT system and the final prosthesis. The prosthesis must be planned under the responsibility of a practitioner, in accordance with the instructions given in the Information for the User. Specific indications for each type of prosthetic component are given below.

straight abutments: suitable for making cemented and screwed prostheses. To make screwed prostheses use the restoration screw with cap SM Q4,5-5;

- 17° abutments**: ideal for cemented prostheses;
- milling abutments**: suitable for making both abutments for prostheses with conometric connection and angled abutments for cemented prostheses, in cases in which 17° abutments cannot be used. Plan the treatment considering that the loading conditions for abutments with high angles put more strain on the implant. The use of milling abutments for making abutments with high angles MUST be compensated by other elements in the restoration treatment;

- Premilled abutments**: suitable for making customized abutments with Medentika® Attachment and maximum working height of 16 mm, using CAD-CAM technique. Suitable for making cemented or screwed prostheses, using automatic milling machines. They are available in one cylindrical section: 11,5 mm (for

angulation to 17°). **NB**: tighten and loosen intact Premilled abutments using ONLY the specific screwdriver Premilled;

- multi abutments**: multifunctional components, suitable for temporary or permanent prostheses, according to the method most suited to the clinical case; available in the **engaging** version and in the **non-engaging** version, free from anti-rotational constraints, to ease insertion even without disparallelism;
- multi abutment castable sleeves**: to combine with multi-abutment to make permanent prostheses with adhesive bonding system, to obtain total passivation of the secondary structures;
- aesthetic bases**: to combine with the use of CAD-CAM systems to make personalised prostheses, they allow the creation of permanent cemented or screw-fitted prostheses with a high level aesthetic quality, guaranteeing a titanium coupling with the implant. To obtain a valid restoration solution do not modify the aesthetic bases;
- temporary aesthetic abutments**: composed of an aesthetic base in titanium and a temporary abutment peek sleeve (that can also be ordered as a spare part). These are suitable for temporary cemented or screw-fitted prostheses of a high aesthetic quality, while guaranteeing a titanium coupling with the implant. The PEEK component must not remain in the oral cavity for more than 180 days;
- temporary abutment peek sleeve**: suitable as replacement for aesthetic base temporary abutments;
- castable abutments**: suitable for making cemented prostheses and ONLY in cases where preformed prosthetic components cannot be used;
- cementable castable abutments**: suitable for making cemented prostheses and ONLY in cases where preformed prosthetic components cannot be used.

NOTE: we recommend using biocompatible casting alloys for all castable components and parts.

During the planning phase of the prostheses, Prodent Italia recommends checking the geometrical relationship between the implant and the prosthesis and to schedule and carry out periodic checks after loading the prosthesis, to verify the stability of the implant, the prosthetic component and the prosthesis, and the correct clinical condition.

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Do not use restoration components to make prostheses other than those specifically indicated for each type. The devices must be used with the sole purpose of completing an implant-restoration treatment made using with endosseous implants belonging to the CONIK-FIT system. The contraindications for endosseous implants of the CONIK-FIT system are listed in the corresponding Information for the User.

4. SIDE EFFECTS AND COMPLICATIONS

After implant-restoration surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. Despite the high success rate in terms of retention of prosthetic implants, the risk of failure may not be completely excluded. The long-term risks that can be associated to implant-restoration include bone reabsorption, failure of integration, oedema, chronic pain and dehiscence. The success of an implant-restoration treatment is strictly connected to how correctly all phases leading to its completion are carried out, including the planning and realisation of the prosthesis. Correct and careful use of the restoration components is a key factor for reducing the risks of implant-restoration treatment failure, in addition to reducing the severity of temporary side effects. Practitioners are asked to report any cases of failure to Prodent Italia S.r.l

5. PACKAGING AND STERILITY

All devices are supplied NON-STERILE, individually packed in a plastic container and inserted in sealed secondary packaging, where the fastening screw is required as accessory to the prosthetic component. It is inserted in the component itself. The label identifying each device is on the secondary packaging. The packing is not suitable for sterilisation; before sterilising the devices, remove them from the packaging and clean them, as indicated in the following paragraph.

6. CLEANING AND STERILISATION

All devices are SINGLE USE. Immediately before use on the patient, they MUST be cleaned and sterilised, using a validated method and under the practitioner's responsibility.(Prodent Italia recomends complying with standard UNI EN ISO 17665 for the development of validation and routine control of the steam sterilisation process, and with standard UNI EN 13060 for determining the testing methods for steam sterilising equipment and for defining the sterilisation cycles to be used). We recommend washing in an ultrasonic bath using appropriate detergents, avoiding using hydrogen peroxide, oxidising acids (such as oxalic, sulphuric and nitric acids) and high concentrations of chlorine, which can damage the devices. Prodent Italia S.r.l. provides further information on the washing and sterilising processes in the *"Important recommendations for the cleaning and sterilisation of devices"* section of the relevant protocol. Prodent Italia S.r.l. has tested the efficacy, on devices of its own manufacture, of the sterilisation process in a saturated steam autoclave at 134°C for 5 minutes.

7. INSERTION AND TIGHTENING INSTRUCTIONS

Tighten all restoration components, except for castable components, exclusively with the specific fastening screw and torque wrench set to 30 Ncm. Tighten castable abutments exclusively with the dedicated fastening screw using the hex screwdriver, without forcing.

NOTE: DO NOT cement the fastening screws inside the abutments; to increase grip we recommend using a dental thread sealant.

NOTE: When using a fastening screw, DO NOT use the torque wrench for tightening components onto the implant analogues.

WARNINGS

- BEFORE** using Restoration Components read carefully the updated Information for the User.
- The name on the device label alone may not be sufficient for identifying its intended use: see Information for the User.
- Restoration Components MUST ONLY be used by qualified dental practitioners who are fully knowledgeable with the theoretical and practical aspects of implantology. Laboratory jobs must be carried out under the instruction and responsibility of the dental practitioner, exclusively by qualified orthodontic technicians.
- CONIK-FIT prosthetic components and the relative fastening screws MUST be used EXCLUSIVELY on implants belonging to the CONIK-FIT implant-restoration system and with the specific prosthetics tool kit.
- When planning the treatment, the dental practitioner is responsible for selecting the implant and the type of prosthetic restoration; the practitioner is likewise responsible for the choice of the type and size of the restoration components and for their correct use.
- ONLY use components of the same prosthetic range as the implant, especially when inserting several implants during the same session, to prevent switches or mistaken implant-screw-abutment couplings.
- Before use, ALWAYS check that the pack, blister and device are undamaged and NEVER use the device if it shows visible signs of damage.
- DO NOT grind or modify the thread and hex screwdriver connections of the fastening screws. This would jeopardise their correct function.
- After the orthodontist has handled the restoration components, before using them on patients ALWAYS check the dimensions and shape of the connection to verify that they have not been altered.
- Take the greatest care during the washing and sterilisation phases, which MUST be performed using a validated method by and under the responsibility of the practitioner.
- Prosthetic components and fastening screws MUST NEVER be reused. Reuse of the device involves a high risk of infection, contamination and implant failure.
- NEVER force the tightening of the devices using the hex screwdriver.
- DO NOT use temporary abutments to make permanent prostheses or prostheses destined to remain in the oral cavity for more than 180 days.
- We recommend that the prosthesis be made in such a way that permanent contact between different metals in the oral cavity are avoided.
- Handle all the devices with the greatest care, to prevent accidental damage that could lead to their incorrect functioning, and during implant-restoration surgery always handle the devices in aseptic conditions and wearing sterile gloves.
- The practitioner is responsible for: diagnosis; treatment planning; the surgical and prosthetic technique adopted; cleaning, sterilisation and handling of the restoration components and relative fastening screws; checking the osteointegration and the stability of the implant and the restoration superstructures; planning and performance of regular check-ups. Mistakes and inaccuracy during these phases may lead to a failure of the implant-restoration treatment.
- The device MUST be handled as hospital waste in the event of disposal. The device should therefore only be disposed of after regular cleaning, in a controlled manner and according to applicable legislation, whenever necessary.
- Prodent Italia S.r.l. ONLY accepts returned devices in their undamaged original packaging and within the contractual terms provided.
- Prodent Italia S.r.l. declines all responsibility for failure to comply with the updates to the Information for the User.
- Prodent Italia S.r.l. furthermore declines all express or implicit liability in relation to direct, indirect or any other type of damage which may be connected to or arise from any evaluation or professional practice error, made during the use of the various products manufactured by Prodent Italia. The use of dental implants, all surgical treatments and prosthetic restorations must always be made under the practitioner's control, who therefore assumes full liability thereof. The user of Prodent Italy's items is solely responsible for evaluating whether a product is suitable or not for each individual patient and clinical case.

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