**IU051-13** – 2021-09-30

**Prodent Italia**

Via Pitagora, 9 – 20016 Pero (MI) Italy – www.prodentitalia.eu

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**ENGLISH**

## CLASS > I INSTRUMENTS **CE** **0546:**

<b>For dental use</b>
<b>INSTRUCTIONS FOR USE</b>
<ol style="list-style-type: none"><li><b>GENERAL INFORMATION</b></li></ol>

All surgical instruments and accessories are made of stainless steel and are invasive, reusable and supplied non-sterile. In order to use the devices correctly, it is necessary to carefully read and apply the Information for the User provided by Prodent Italia. This information includes these Instructions for Use and the relevant clinical/surgical protocol for each device, as indicated in the INDICATIONS paragraph below. The Information for the User can be found on the company website or by contacting Prodent Italia S.r.l. directly.

- INDICATIONS**

These instruments are intended for use when preparing the surgical site and inserting Prodent Italia endosseous devices in the implant site. The devices are invasive and MUST be cleaned and sterilised immediately before use. The devices can either be specific to a certain implant range, or to a certain implant-restoration connection or common to several ranges or connections. Some specific indications for use are provided below; the appropriate protocol is specified if the device is NOT available for all implant types. If the device is available for all implant types, refer to the specific protocol for the implant being inserted: CL007 for PROWIDE implants, CL008 for PRIME and TWINNER implants, CL010 for CONOMET implants and CL011 for PROSHAPE implants. The devices identified as PROFAST are intended for use with PROFAST implants ONLY, refer to clinical protocol CL012.

- **Contra-angle hex screwdriver**: Available with two different lengths, ONLY to be used connected to the contra-angle to mechanically tighten and loosen the fastening screws of prosthetic components with a hexagonal recess at 30 Ncm. **N.B.:** The devices have 2 laser-etched notches on the shaft between the tip and the handle, to distinguish them from the previous version screwdrivers, which are not marked. **N.B.:** Do NOT use unmarked screwdrivers to tighten devices with TS connections (CONOMET and PROSHAPE implants). **N.B.:** For those prosthetic components that require it, subsequently perform final tightening manually with the hex bit for the torque wrench.

- **Driver**: Available as a single variant for each restoration range and for each connection; to be used for picking up the implants, positioning them in the receptor site and starting treatment. Insert the driver connected to the digital wrench in the implant connection. Check that the driver has completely engaged the implant connection; if necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring marked on the driver as reference). Move the implant to the receptor site, position it so that it is firm and stable and disconnect the digital wrench. Connect the driver to the Torque wrench and tighten the implant until it is correctly positioned. Follow the torque indications provided in the protocols and Instructions for Use for the implants to be used. The drivers are fitted with O-rings and a dedicated square connector for connection with the digital wrench and the torque wrench. To identify the size variants, refer to the specific protocol for the implant being inserted.

- **Contra-angle Driver**: Available in different variants for different restoration ranges; to be used connected to the contra-angle for picking up the implants, positioning them in the receptor site and inserting them without EVER exceeding 25 rpm and 45 Ncm; Insert the driver connected to the contra-angle in the implant connection. Check that the driver has completely engaged the

implant connection; if necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring marked on the driver as reference). Take the implant to the receptor site and tighten the implant until it is correctly positioned. To identify the size variants, refer to the specific protocol for the implant being inserted.

- **Extension**: Available as a single variant; to be used when it is deemed necessary to increase the length of the connection between the digital wrench or the torque wrench or connection for the contra-angle and dedicated devices. It is fitted with an O-ring for connection with the above devices.

- **Bur extension**: Available as a single variant (Ref. Clinical Protocols); to be used when it is deemed necessary to extend the connection between the bur and the contra-angle. Do not use the device beyond 45 Ncm.

- **Contra-angle connection**: Available as a single variant; to be used for mechanical use of the screw taps and for mechanical implant insertion. NEVER exceed 25 rpm and 45 Ncm during use.

- CONTRAINDICATIONS**

All instruments must be used exclusively according to the indications provided and to prepare the site and insert the devices they are to be used for. The contraindications concerning Prodent Italia implantable devices are listed in the corresponding Information for the User.

- SIDE EFFECTS AND COMPLICATIONS**

After oral surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. The long-term risks that can be associated with implant-restoration or guided bone regeneration procedures include bone resorption, unsuccessful integration, oedema, chronic pain and dehiscence. The success of implant-restoration or guided bone regeneration treatment is closely related to the correctness of all the phases involved, including the preparation of the receptor site. Correct and careful use of the instruments is a key factor for reducing the risks of implant-restoration or guided bone regeneration treatment failure, in addition to reducing the severity of temporary side effects. Practitioners are asked to report any cases of failure to Prodent Italia S.r.l

- PACKAGING AND STERILITY**

All devices are provided NON-STERILE and individually packed; depending on the nature of each one, they may be supplied in a plastic container inside sealed secondary packaging, or directly inside heat-sealed transparent plastic bags or directly inside a sealed box. The label identifying each device is placed, according to the cases, on the secondary packaging, box or plastic bag. No packaging is suitable for sterilisation. Before sterilising the devices, remove them from the packaging and clean them, as indicated in the following paragraph.

- CLEANING AND STERILISATION**

All devices are reusable. Immediately before use, they MUST be **cleaned and sterilised, using a validated method and under the practitioner’s responsibility** (Prodent Italia recommends complying with standard UNI EN ISO 17665 for the development of validation and routine control of the steam sterilisation process, and with standard UNI EN 13060 for determining the testing methods for steam sterilising equipment and for defining the sterilisation cycles to be used). Immediately after use, the devices MUST be cleaned to remove all types of residue, using plastic brushes with stiff, non-metal bristles. We recommend washing in an ultrasonic bath using appropriate detergents, avoiding using hydrogen peroxide, oxidising acids (such as oxalic, sulphuric and nitric acids) and high concentrations of chlorine, which can damage the devices. Prodent Italia S.r.l. provides **further information on the washing and sterilising processes in the “Important recommendations for the cleaning and sterilisation of devices”** section of the relevant protocol. Prodent Italia S.r.l. has tested the efficacy, on devices of its own manufacture, of the sterilisation process in a saturated steam autoclave at 134° C for 5 minutes.

**⚠ WARNINGS:**

- ONLY use the devices after reading the most recent version of the Information for the User.
- The name on the device label alone may not be sufficient to identify its intended use: refer to the Information for the User.
- All devices MUST ONLY be used-by qualified dental professionals who are fully knowledgeable with the theoretical and practical aspects of surgical procedures and under their own responsibility.
- All devices must only be used according to the indications provided and with the devices they are intended for. When available as different sized variants, identify using markings and colour-coding.
- Use the devices in question with original Prodent Italia S.r.l. devices only.
- NEVER use devices with damaged markings or colour-coding.
- Before first use ALWAYS check that the device pack is intact. NEVER use devices if the pack is visibly damaged.
- ALWAYS check the device is undamaged before each use. NEVER use devices if they are visibly damaged. If in doubt as to the integrity of the devices, request a check by Prodent Italia.
- Take the greatest care during the washing and sterilisation phases, which MUST be performed using a validated method by and under the responsibility of the practitioner.
- NEVER use devices fitted with an O-ring, if it is missing or worn.
- Before use, it is the practitioner’s responsibility to make sure that each device is fit for use and suitable for the intended use and for the envisaged purposes, also according to the Information for the User.
- When connecting two or more devices, ALWAYS check that they are correctly and fully connected before and during use as specified in the indications for use.
- All devices should be handled with special care, to avoid accidental damage that could jeopardise the results of treatment.
- During use on patients, handle all devices in aseptic conditions and wearing sterile gloves.
- NEVER apply excessive force when tightening.
- The practitioner is responsible for: the planning of the implant-restoration treatment; the surgical technique adopted; the choice and size of device to be implanted; the choice of rehabilitation technique; the choice and use of restoration instruments; the creation of the restoration and the planning and performance of regular check-ups. Mistakes and inaccuracy during these phases may lead to the failure of the implant-restoration treatment.

- The device MUST be handled as hospital waste in the event of disposal. The device should therefore only be disposed of after proper cleaning, in a controlled manner and according to current applicable legislation, whenever necessary.
- Prodent Italia S.r.l. ONLY accepts returned devices in their undamaged original packaging and within the contractual terms provided.
- Prodent Italia S.r.l. declines all responsibility for failure to comply with the updates to the Information for the User.
- Prodent Italia S.r.l. furthermore declines all express or implicit liability in relation to direct, indirect or any other type of damage which may be connected to or arise from any evaluation or professional practice error, made during the use of the various products manufactured by Prodent Italia. The use of dental implants, all surgical treatments and prosthetic restorations must always take place under the supervision of the practitioner, who therefore assumes full liability. Users of Prodent Italia products have sole responsibility for evaluating whether a product is suitable for each individual patient and clinical case.

**KEY for the symbols present on the label:**

<b>LOT</b>	Batch code.
<b>REF</b>	Catalogue number.
<b>⚠</b>	Caution! The device is associated with warnings not indicated on the label and that are described in the Instructions for Use.
<b>i</b>	Consult the Instructions for Use inside the pack.
<b>⊘</b>	DO NOT use if the package is damaged.
<b>⚙</b>	Manufacturer pursuant to Directive (EU) 93/42/EEC as amended.

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- Utilizați dispozitivele din obiect numai cu alte dispozitive originale Prodent Italia S.r.l.