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### INSTRUCTIONS FOR USE

#### 1. GENERAL INFORMATION

The healing screws and caps are made of grade 5 titanium and are invasive, reusable devices, that are supplied non-sterile. In order to use the devices correctly, it is necessary to carefully read and apply the Information for the User provided for each of them by Prodent Italia. Such information consists of these Instructions for Use and the clinical/surgical protocol relative to the implants to be inserted, i.e. protocols CL008 for PRIME and TWINNER implants, CL010 for CONOMET implants, CL011 for PROSHAPE implants, CL007 for PROWIDE implants and CL012 for PROFAST implants. This information can be found online in the company website or by contacting Prodent Italia S.r.l. directly. The protocols indicate the types and sizes of healing screws and caps that are available to meet the various anatomical and healing requirements. The healing screws and caps are identified by laser markings and/or colours, as indicated in the protocols.

#### 2. INDICATIONS

The healing screws and caps are indicated for transmucosal healing of the tissue and for conditioning the gums, before fitting the prosthesis to the implant. In the case of PROFAST range, suitable for immediate loading on more implants, the healing screw is also useful to protect the transmucosal connection during the surgery and in the phases prior to the restoration. All devices can be screwed: tighten the **healing screws** in the implant connection; tighten the **FAST healing cap** to the FAST base already screwed into the implant connection; tighten the **PROFAST healing cap** directly to PROFAST implant. Tighten all devices **EXCLUSIVELY** using a hex screwdriver. To achieve the correct healing and gum conditioning, only use healing screws of GREATER height than the height of the mucous membrane to be conditioned. To achieve the correct gum conditioning, use healing screws with the same transmucosal section as the prosthesis destined to fit on the implant. Cylindrical healing screws are suitable in case of restoration with implant components having cylindrical transmucosal section and in case there is so little space between the implants and/or the residual dentition that tapered healing screws cannot be used. Prodent Italia recommends assessing the primary implant stability before inserting the healing screws or caps. It also recommends advising the patient not to chew on the healing screws or caps, and to schedule periodic checks to evaluate the stability of the implant and its correct clinical condition.

#### 3. CONTRAINDICATIONS

All healing screws and caps must be used exclusively in accordance with the instructions and only with the endosseous implants they are dedicated to. The contraindications for implants are listed in the corresponding Information for the User.

#### 4. SIDE EFFECTS AND COMPLICATIONS

After dental implant surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. The long-term risks that can be associated to implant-restoration include bone reabsorption, failure of integration, oedema, chronic pain and dehiscence. The success of the implant treatment is strictly connected to how correctly all phases leading to its completion are carried out, including the healing of the bone tissue and the gum conditioning. Correct and careful use of the restoration components is a key factor for reducing the risks of implant-restoration treatment failure, in addition to reducing the severity of temporary side effects. Practitioners are asked to report any cases of failure to Prodent Italia S.r.l.

#### 5. PACKAGING AND STERILITY

All devices are supplied NON-STERILE, individually packed in a plastic container and inserted in sealed secondary packaging. The label identifying each device is on the secondary packaging. The packing is not suitable for sterilisation; Before sterilising the devices, remove them from the packaging and clean them, as indicated in the following paragraph.

#### 6. CLEANING AND STERILISATION

All devices are reusable. Immediately before use, they **MUST be cleaned and sterilised, using a validated method and under the practitioner's responsibility** (Prodent Italia recommends compliance with standard UNI EN ISO 17665 for the development of validation and routine control of the steam sterilisation process, and standard UNI EN 13060 for determining the testing methods for steam sterilising equipment and for defining the sterilisation cycles to be used). Immediately after use, the devices **MUST** be cleaned to remove all types of residue, using plastic brushes with stiff, non-metal bristles. We recommend washing in an ultrasonic bath using appropriate detergents, avoiding using hydrogen peroxide, oxidising acids (such as oxalic, sulphuric and nitric acids) and high concentrations of chlorine, which can damage the devices. Prodent Italia S.r.l. provides detailed information on the washing and sterilising processes in the "Important recommendations for the cleaning and sterilisation of devices" section of the relevant protocol. Prodent Italia S.r.l. has tested the efficacy, on devices of its own manufacture, of the sterilisation process in a saturated steam autoclave at 134°C for 5 minutes.

#### WARNINGS

- ONLY use the healing screws and caps after carefully reading the updated information for the User.
- The name on the device label alone may not be sufficient for identifying its intended use: see Information for the User.
- The healing screws and caps **MUST ONLY** be used by qualified dental practitioners who are fully knowledgeable with the theoretical and practical aspects of implantology.
- The healing screws and caps **MUST** be used **ONLY** on implants belonging to the same implant-restoration system and with specific instruments.
- **ALWAYS** check that the packaging is undamaged before using the device the first time. **DO NOT** use healing screws and caps if their packaging shows signs of damage.
- **ALWAYS** check the device is undamaged before each use. **DO NOT** use healing screws and caps if they show signs of damage. If in doubt as to the integrity of the devices, request a check by Prodent Italia.
- Take the greatest care during the washing and sterilisation phases, which **MUST** be performed using a validated method by and under the responsibility of the practitioner.
- Before use, it is the practitioner's responsibility to make sure that the device is fit for use and suitable for the intended use and for the envisaged purposes, also according to the Information for the User.
- **ONLY** use healing screws and caps belonging to the same range as the implant, paying attention to the markings and/or colour, especially in cases where several implants are inserted during the same session, to prevent switches or mistaken implant-screw coupling (See protocols CL007, CL008, CL010, CL011).
- **DO NOT** use healing screws and caps if the markings or colour are deteriorated.
- All devices should be handled with special care, to avoid accidental damage that could jeopardise the results of treatment.
- Handle healing screws and caps in aseptic conditions and with sterile gloves.
- **NEVER** force the tightening of the devices using the hex screwdriver.
- The practitioner is responsible for planning the implant-restoration treatment, the surgical technique adopted, the choice and size of the device to be implanted, the choice and use of the healing screws and the caps, the choice of rehabilitation technique as well as scheduling and performance of routine checks. Mistakes and inaccuracy during these phases may lead to a failure of the implant-restoration treatment.
- The device **MUST** be handled as hospital waste in the event of disposal. The device should therefore only be disposed of after regular cleaning, in a controlled manner and according to applicable legislation, whenever necessary.
- Prodent Italia S.r.l. **ONLY** accepts returned devices in their undamaged original packaging and within the contractual terms provided.
- Prodent Italia S.r.l. declines all responsibility for failure to comply with the updates to the Information for the User.
- Prodent Italia S.r.l. furthermore declines all express or implicit liability in relation to direct, indirect or any other type of damage which may be connected to or arise from any evaluation or professional practice error, made during the use of the various products manufactured by Prodent Italia. The use of dental implants, all surgical treatments and prosthetic restorations must always be made under the practitioner's control, who therefore assumes full liability thereof. The user of Prodent Italia's items is solely responsible for evaluating whether a product is suitable or not for each individual patient and clinical case.

#### KEY for the symbols present on the label:

Batch code.

Catalogue number.

Caution! The device is associated with warnings not indicated on the label and that are described in the Instructions for Use.

Consult the Instructions for Use inside the pack

DO NOT use if the package is damaged

Manufacturer pursuant to 93/42/EEC and subsequent amendments