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## ENGLISH

### PREFORMED AND SEMI-CASTABLE PROSTHETIC COMPONENTS 0546

For dental use

#### INSTRUCTIONS FOR USE

##### 1. GENERAL INFORMATION

The prostheses are SINGLE-USE, they are supplied non-sterile and are available in various types and sizes to create prostheses that meet different anatomical and functional requirements. The prostheses include preformed and semi-castable prosthetic components. All devices are invasive. The devices to be connected to the implant using fastening screws are supplied together with the dedicated screw. The table lists the types of prosthetic components, the materials they are made of and the dedicated accessories they are supplied with.

Restoration components and Screws	Dedicated accessories (provided in pack)	Material
<b>PREFORMED</b>		
Ball attachments	Insertion support in acetal resin	Grade 5 titanium
FAST 0° bases		
FAST angled bases		
Straight abutments (0°)		
Angled abutments (17° - 25°)		
Milling abutments		
Multi abutment		
Aesthetic bases		
Link bases		
Connect bases		
FAST temporary abutments	Fastening screw	Grade 5 titanium base PEEK bush
PROFAST temporary abutments		
Premilled abutments		
Torx Screws for angulated holes		
Temporary aesthetic base abutments	Fastening screw	Grade 5 titanium base PEEK bush
Temporary abutment peek sleeve	/	PEEK
Temporary abutments	Fastening screw	Acetal resin
<b>SEMI-CASTABLE</b>		
Rod abutments	Fastening screw	Grade 5 titanium base Polycarbonate castable cap
Restoration screws with cap	/	Grade 5 titanium screw, Polycarbonate cap

All the devices, to be used correctly, require practitioners to carefully read and apply the Information for the User provided by Prodent Italia. This information is contained in these Instructions for Use and in the specific Clinical protocol for each type of implant, which is indicated in the INDICATIONS section. All information for Users can be found on the company website or by contacting Prodent Italia S.r.l. directly.

All Fastening Screws for prosthetic components are made of Grade 5 titanium, intended for SINGLE-USE and supplied non-sterile. The specific fastening screws for the different prosthetic components are listed in the corresponding Clinical Protocols. The clinical protocols also indicate the types and sizes of the prosthetic components and corresponding fastening screws to meet the various anatomical and healing requirements. The prosthetic range of each component is identified by markings and/or colour coding as indicated in the Protocol, except for the prosthetic components for PROWIDE implants, that are NOT marked or colour-coded.

##### 2. INDICATIONS

The restoration components are interconnection devices for the endosseous implants and the prosthesis. The prosthesis may be a single crown, a bridge, a partial prosthesis or an overdenture, according to the treatment the dental practitioner has planned, under his own responsibility. Specific indications are given below for each type of prosthetic component, which may not be available for all implant types. In order to correctly identify the types of implant that each single device is indicated for, users should consult the specific Clinical Protocols for each type of implant: CL007 for PROWIDE implants, CL008 for PRIME and TWINNER implants, CL010 for CONOMET implants and CL011 for PROSHAPE implants, CL012 for PROFASST implants. **N.B.:** PROWIDE implants have a CONIK-FIT connection, PRIME and TWINNER implants have a PERFORM SM connection; CONOMET and PROSHAPE implants have a TS connection. PROFASST implants have a single transmucosal connection for all implant sizes.

- **Straight abutments (0°):** suitable for cemented prostheses; in the case of abutments with a coronal height of 3 mm ONLY, they can be used in combination with restoration screws with caps to achieve screw-retained prostheses.

- **N.B.** DO NOT grind abutments with a coronal height of 3 mm if they are to be used with restoration screws with caps in order to achieve a screw-retained prosthesis.

- **Angled abutments (17° - 25°):** ideal for cemented prostheses;

- **Milling abutments:** Indicated for cemented and conometric prostheses; DO NOT use for angled abutments when the divergence between pillars is more than 20°.

- **Premilled abutments:** suitable for making customized abutments with Medentika® Attachment and maximum working height of 16 mm, using CAD-CAM technique. Suitable for making cemented or screwed prostheses, using automatic milling machines. They are available in two cylindrical sections: 11,5 mm (for angulation to 17°) and 15,8 mm (for angulation to 25°). **N.B.:** tighten and loosen intact Premilled abutments using ONLY the specific screwdriver Premilled.

- **Multi abutments:** multifunctional components, suitable for temporary or permanent prostheses, according to the method most suited to the clinical case; available in the **engaging** version and in the **non-engaging** version, free from anti rotational constraints, to ease insertion even in the presence of disparallelism.

- **Aesthetic bases and connect bases:** to combine with the use of CAD-CAM systems to make personalised prostheses with a bonding technique, they allow the creation of permanent cemented or screw-fitted prostheses with a high-level aesthetic quality, guaranteeing a titanium coupling with the implant. To obtain a valid restoration solution do not modify the aesthetic bases. Connect bases are available with an implant-restoration connection in the **engaging** version and in the **non-engaging** version, free from anti rotational constraints, to ease insertion even in the presence of disparallelism; each version is available in different variants with different transmucosal heights (see specific Clinical Protocol).

**N.B.:** the transmucosal section that can be obtained using the Connect bases is NOT the same as that of the healing screws; in order to obtain a valid prosthetic solution with a good cosmetic result, before fitting the final restoration, prime the gum with a temporary restoration using a Connect base.

- **Link bases:** for the creation of customised prostheses using a bonding technique; they can be used to achieve cement- or screw-retained permanent restorations with an excellent aesthetic result, whilst guaranteeing a titanium coupling with the implant. The coronal portion can be cut as indicated in the Clinical Protocols to achieve a height suited to the specific clinical case. Link TS bases are available as an **engaging**

version and a **non-engaging** version, that is free from antirrotational constraints, to ease insertion even in the presence of disparallelism. Link SM bases are available as an engaging version only. Link bases with different transmucosal heights are available for each type of connection (see the specific Clinical Protocols). Link bases can be fitted either with a special Scanmarker for Link Bases to be retained using a dedicated screw in order to take an impression digitally or with an EasyTransfer device in order to take an impression on a single tooth using a closed tray technique with a non-perforated impression tray.

**N.B.:** before creating a permanent restoration with a Link TS base, the tissues must be conditioned using the dedicated TS Link healing screws.

**N.B.:** to use Scanmarkers for Link bases correctly, fit the Scanmarker on to the Link base without a screw, then fit the base together with the Scanmarker on to the implant or lab analogue and lastly secure the whole assembly using the dedicated Scanmarker screw.

**N.B.:** Do NOT grind Link bases before using them to take impressions with EasyTransfer.

- **Temporary aesthetic abutments:** composed of an aesthetic base in titanium and a temporary abutment peek sleeve (that can also be ordered as a spare part). These are suitable for temporary cemented or screw-fitted prostheses of a high aesthetic quality, while guaranteeing a titanium coupling with the implant. The Peek component must not remain in the oral cavity for more than 180 days.

- **Temporary abutment peek sleeves:** usable as replacements for aesthetic base temporary abutments.

- **Temporary abutments:** suitable for temporary cemented prostheses, must not be left in oral cavity more than 30 days.

- **FAST and PROFASST temporary abutments:** suitable for making temporary prostheses respectively on FAST bases and on PROFASST implants; **N.B.:** Do NOT cut the abutments under the first coronal ring from the cone.

- **FAST Bases:** suitable for making screwed total prostheses; may be used screwed to the implant ONLY as support structure for the specific FAST temporary abutments and FAST castable abutments; DO NOT use grind; DO NOT use for prostheses on a single implant; DO NOT use the FAST 0° bases if the divergence between the implants exceeds 6°.

- **Ball attachments:** suitable for holding overdenture removable prostheses resting on gums, coupled with specific Normal size caps and containers; DO NOT use for prostheses on a single implant; DO NOT use if implant divergence angle exceeds 10°; DO NOT grind. These are supplied with a specific insertion support, which can be used as impression transfer, as described in the specific section of the Clinical Protocol.

- **Rod abutments:** suitable for making removable, partial or overdentures prostheses;

- **Restoration screws with cap:** use EXCLUSIVELY on straight abutments in titanium with crown height of 3 that have not been ground for making screwed prostheses; see Protocol CL008.

- **Torx screws for angulated holes:** intended for securing restorations for which the hole is at an angle to the axis of the implant, using the dedicated screwdriver. The screws have a Ø2.3 head with a torx recess. **N.B.:** torx screws can only be used with customised restorations produced considering the size of the screw or with Link bases. DO NOT use with other types of preformed or semi-castable prosthetic parts. **N.B.:** tighten torx screws using a dynamometric key set to 20 Ncm and, in any case, NEVER to more than 25 Ncm.

**N.B.:** We recommend using biocompatible casting alloys for all castable parts.

**N.B.:** The devices belonging to FAST and PROFASST ranges are also suitable for immediate loading on temporary prostheses; the dental practitioner is responsible for adequately selecting the patient, the surgical procedure and the devices applied, and for checking the primary stability of the implant prior to applying load.

All prosthetic components, except for the FAST bases, the Ball Attachments, the Aesthetic Bases and the Connect Bases, may be modified by qualified dental staff, to make the required prosthesis. During the planning phase of the prostheses, Prodent Italia recommends checking the geometrical relationship between the implant and the prosthesis and to schedule and carry out periodic checks after loading the prosthesis, to verify the stability of the implant, the prosthetic component and the prosthesis, and the correct clinical condition.

##### 3. CONTRAINDICATIONS

Do not use restoration components to make prostheses other than those specifically indicated for each type. Prosthetic parts with non-engaging connections are contraindicated for single-implant restorations. The devices must be used for the sole purpose of completing an implant-restoration treatment with endosseous components provided with the corresponding connection. The contraindications for the endosseous implants are listed in the corresponding Information for the User.

##### 4. SIDE EFFECTS AND COMPLICATIONS

After implant-restoration surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. Despite the high success rate in terms of retention of prosthetic implants, the risk of failure may not be completely excluded. The long-term risks that can be associated to implant procedures include bone resorption, failure of integration, oedema, chronic pain and dehiscence. The success of an implant-restoration treatment is strictly connected to how correctly all phases leading to its completion are carried out, including the planning and realisation of the prosthesis. Correct and careful use of the restoration components is a key factor for reducing the risks of implant-restoration treatment failure, in addition to reducing the severity of temporary side effects. Practitioners are asked to report any cases of failure to Prodent Italia S.r.l.

##### 5. PACKAGING AND STERILITY

All devices are supplied NON-STERILE, individually packed in a plastic container and inserted in sealed secondary packaging. Where the fastening screw is required as accessory to the prosthetic component, it is inserted in the component itself. The label identifying each device is on the secondary packaging. The packing is not suitable for sterilisation; Before sterilising the devices, remove them from the packaging and clean them, as indicated in the following paragraph.

##### 6. CLEANING AND STERILISATION

All devices are SINGLE-USE. Immediately before use, they MUST be cleaned and sterilised using a validated method and under the practitioner's responsibility. (Prodent Italia recommends compliance with standard UNI EN ISO 17665 for the development of validation and routine control of the steam sterilisation process, and with standard UNI EN 13060 for determining the testing methods for steam sterilising equipment and for defining the sterilisation cycles to be used). We recommend washing in an ultrasonic bath using appropriate detergents, avoiding using hydrogen peroxide, oxidising acids (such as oxalic, sulphuric and nitric acids) and high concentrations of chlorine, which can damage the devices. Prodent Italia S.r.l. provides further information on the washing and sterilising processes in the "Important recommendations for the cleaning and sterilisation of devices" section of the relevant protocol. Prodent Italia S.r.l. has tested the efficacy, on devices of its own manufacture, of the sterilisation process in a saturated steam autoclave at 134°C for 5 minutes.

##### 7. INSERTION AND TIGHTENING INSTRUCTIONS

The components directly attachable to the implant (ball attachments and FAST 0° bases) are supplied with an insertion support that allows them to be positioned in the implant connection before the final tightening. Tighten straight abutments, angled abutments, milling abutments, multi abutments, aesthetic bases, temporary aesthetic abutments, FAST angled bases ONLY with the specific fastening screws and to 30 Ncm using a torque settable instrument. If Link bases are used with torx screws, they should be tightened at 20 Ncm and, in any case, never at more than 25 Ncm. Tighten

FAST and PROFASST temporary abutments only with the specific fastening screws to 20 Ncm, using a torque settable instrument. Tighten acetal resin temporary abutments, the rod abutments, and the restoration screws with caps, ONLY using the dedicated fastening screws and a hex screwdriver, without forcing. Tighten the ball attachments and FAST 0° bases using hex wrench CH.2.6 with the torque gauge set to 30 Ncm.

**N.B.:** when tightening prosthetic components with a TS connection (for CONOMET and PROSHAPE implants) ONLY use screwdrivers with a double laser-etched ring on the shaft, do not use screwdrivers without a double ring marking.

**N.B.:** DO NOT cement the fastening screws inside the abutments; to increase grip we recommend using a dental thread sealant.

**N.B.:** When using a fastening screw, DO NOT use the torque wrench for tightening components onto the implant analogues.

##### ⚠️ WARNINGS

- BEFORE using Restoration Components read carefully the updated Information for the User.

- The name on the device label alone may not be sufficient for identifying its intended use; see Information for the User.

- Restoration Components MUST ONLY be used by qualified dental practitioners who are fully knowledgeable with the theoretical and practical aspects of implantology. Laboratory jobs must be carried out under the instruction and responsibility of the dental practitioner, exclusively by qualified dental technicians.

- Pay careful attention to the implant-restoration connection; use prosthetic components on implants with the same connection and using the appropriate restoration instrument kit ONLY; PROWIDE implants have a single CONIK-FIT connection for all sizes; PRIME and TWINNER implants have PERFORM SM connections with different dimensions depending on the implants diameter; the connection diameter defines the restoration range (see CL008); CONOMET and PROSHAPE implants have a single TS connections for all implants sizes; PROFASST implants have a single transmucosal connection for all sizes. Pay careful attention to the markings and/or colours as specified in the protocols, especially in case of insertion of more implants in the same context, to avoid inversions or wrong implant-screw-abutment coupling.

- When planning the treatment, the dental practitioner is responsible for selecting the implant and the type of prosthetic restoration. The practitioner is likewise responsible for the choice of the type and size of the restoration components and for their application.

- Before use ALWAYS check that the pack, blister and device are undamaged and NEVER use the device with visible signs of damage.

- DO NOT grind or modify the thread and hex screwdriver connections of the fastening screws. This would jeopardise their correct function.

- NEVER use torx screws for angulated holes on parts other than those specified in the Indications section.

- After the dental technician has handled the prosthetic components, before using them on patients ALWAYS check the dimensions and shape of the connection to verify that they have not been altered.

- NEVER use prosthetic parts with a non-engaging connection for single-implant restorations.

- Take the greatest care during the washing and sterilisation phases, which MUST be performed using a validated method by and under the responsibility of the practitioner.

- Prosthetic components and fastening screws MUST NEVER BE REUSED. Reuse of the device involves a high risk of infection, contamination and implant failure.

- NEVER force the tightening of the devices using the hex screwdriver.

- NEVER use temporary abutments as castable abutments. NEVER use temporary abutments for making permanent prostheses or prostheses destined to remain in the oral cavity for more than the time indicated for each type, except for that methodologies specified in Instructions for Use.

- When joining together two or more devices, in accordance with the instructions for use, ALWAYS check that they are correctly and completely connected before and during use.

- NEVER use devices with markings or colour-coding if the marking or colour coding is damaged.

- We recommend that the prosthesis be made in such a way that permanent contact between different metals in the oral cavity is avoided.

- Handle all devices with special care, to prevent accidental damage that could lead to their incorrect functioning, and during surgery always handle devices in aseptic conditions while wearing sterile gloves.

- The practitioner is responsible for: diagnosis; treatment planning; the surgical and prosthetic technique adopted; the cleaning, sterilisation and handling of the restoration components and relative fastening screws; checking the osteointegration and the stability of the implant and the restoration superstructures; the planning and performance of regular check-ups. Mistakes and inaccuracy during these phases may lead to a failure of the implant-restoration treatment.

- The device MUST be handled as hospital waste in the event of disposal. The device should therefore only be disposed of after regular cleaning, in a controlled manner and according to applicable legislation, whenever necessary.

- Prodent Italia S.r.l. ONLY accepts returned devices in their undamaged original packaging and within the contractual terms provided.

- Prodent Italia S.r.l. declines all responsibility for failure to comply with the updates to the Information for the User.

- Prodent Italia S.r.l. furthermore declines all express or implicit liability in relation to direct, indirect or any other type of damage which may be connected to or arise from any evaluation or professional practice error, made during the use of the various products manufactured by Prodent Italia. The use of dental implants, all surgical treatments and prosthetic restorations must always be made under the practitioner's control, who therefore assumes full liability thereof. The user of Prodent Italia's items is solely responsible for evaluating whether a product is suitable or not for each individual patient and clinical case.

##### KEY for the symbols present on the label:

	DO NOT reuse! The device is SINGLE-USE and intended to be used once and on one patient only. Reuse could cause serious risks of cross-infection and/or loss of functionality.
	Batch code
	Catalogue number
	Caution! The device is associated with warnings not indicated on the label and that are described in the Instructions for Use.
	Consult the Instructions for Use inside the pack.
	DO NOT use if the package is damaged
	Manufacturer pursuant to 93/42/EEC and subsequent amendments