COMPONENTI ENDOSSEI PROSHAPE C E 0546

Ad uso odontoiatrico

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For dental use

INSTRUCTIONS FOR USE**1. GENERAL INFORMATION**

The PROSHAPE range implants are endosseous dental implants in pure medical-grade titanium available in the following different types: PROSHAPE INTEGRAL and PROSHAPE HYBRID. PROSHAPE implants have a TS implant-restoration connection. Each implant is provided with a dedicated cover screw. All the devices contained in the pack are **STERILE** and **SINGLE-USE**. The endosseous surface varies according to the different type of implant: the HYBRID implants feature a machined section of the neck, while the remaining surface is MPS-treated; the INTEGRAL implants feature an entirely MPS-treated bone contact surface. All implants are available in a number of different diameters and heights, to meet diverse anatomical requirements and, in order to be used correctly, practitioners must carefully read and apply the information for the User provided by Prodent Italia. This information consists in the Instructions for Use and Clinical Protocol CL011. The information for the User can also be found in the company website or by contacting Prodent Italia S.r.l. directly.

2. INDICATIONS

The different types of PROSHAPE implants, considered as a whole, allow us to treat all cases of aesthetic and functional rehabilitation using endosseous dental implants. They can be used in both the upper and lower arches in patients with partial or total edentulism. Restoration may involve single crowns, bridges, partial or full-arch dentures. HYBRID implants are primarily indicated when the practitioner, having evaluated the specific clinical case, decides that part of the implant neck should protrude from the bone crest. Prodent Italia recommends that you evaluate the bone quality during surgery planning, the size of the implant to be used depending on the position in the oral cavity, the implant primary stability at the end of placement, the implant-to prostheses geometric ratio during restoration planning, and that you schedule and conduct regular check-ups. The general instructions to be followed by the patient after surgery are given on the back of the Dental Implant Passport, found inside each implant pack.

3. CONTRAINDICATIONS

Endosseous implants are contraindicated in patients who, in general: are in a poor state of general health; have a systemic illness or cancer; have severe immune system, cardiac, neurological or endocrine disorders; have chronic or acute infectious diseases; have serious vascular problems; have chronic jaw osteitis; have inadequate and/or poor quality bone tissue; suffer from bruxism; have poor oral hygiene; are addicted to nicotine; have an unsuitable mental profile; are drug users; are alcohol abusers; are pregnant; are allergic to titanium or similar metals; are frequently exposed to radiation. Any adverse periodontal problems must be resolved before placing implants.

4. SIDE EFFECTS AND COMPLICATIONS

After dental implant surgery, patients may experience temporary side effects such as pain, swelling, haematomas, speech problems, and inflammations of the soft tissues. Despite the high success rate for implant procedures and their stability, the risk of failure cannot be completely ruled out. The long-term risks that can be associated to implant procedures include bone reabsorption, failure of integration, oedema, chronic pain and dehiscence. Implant treatment success is closely related to the correctness of the medical diagnosis, treatment planning and the surgical and restoration techniques adopted. Practitioners are asked to report any instance of failure to Prodent Italia S.r.l.

5. PACKAGING AND STERILITY

Each implant is supplied **STERILE** inside a plastic container. The implant is fitted in a pure medical-grade titanium ring. In addition, the apex of the implants rests on a pure medical-grade titanium disk, in turn supported by a plastic spacer. The cover screw is fitted into the cap which closes the implant housing. The plastic container is supplied inside a sealed plastic blister to preserve sterility. The blister comes inside a cardboard box with closing seals for optimal storage. The implant label is

placed on the box. The box also contains the Instructions for Use, the patient's Dental Implant Passport and the extra peel-off label with the device's identification details to be attached to the patient's clinical records. The implant housing and its extraction instruments make it possible to avoid contact between the surface of the implant and other components and surfaces other than titanium before placement in the site.

6. STORING AND HANDLING THE PACK

The sterility of each implant until the expiry date indicated on the label is assured by storing it inside its original pack, sealed and undamaged, in a dry place at room temperature. Before using the implant, always check that the pack is undamaged and has no visible signs of damage that could compromise its sterility. The pack must not be opened until the implant is to be used. The blister must be opened and the implant taken out in conditions of asepsis.

7. INDICATIONS FOR IMPLANT PLACEMENT

After preparing the implant site (see Protocol CL011), open the blister and remove the device plastic container. Open the plastic container by turning the implant housing cap anti-clockwise and pulling upwards to remove it. **N.B.** The cover screw is fitted into the cap of the implant housing.

7.1 MANUAL IMPLANT PLACEMENT

Pick the implant up using a driver connected to the Digital Wrench. Insert the driver in the implant connection. Check that the driver has completely engaged the implant connection. If necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring marked on the driver as reference). Move the implant to the receptor site, position it so that it is firm and stable and disconnect the digital wrench. Connect the driver to the Ratchet and tighten the implant until it reaches its correct final position; NEVER exceed 25 Ncm. When insertion is complete, disconnect the ratchet and driver from the implant.

7.2 MECHANICAL IMPLANT PLACEMENT

Pick the implant up using a driver for contra-angle connected to the contra-angle. Insert the driver in the implant connection. Check that the Contra-angle Driver has completely engaged the implant connection. If necessary, turn the driver to ease insertion. Remove the implant from the container and visually check that there is no gap between the driver and the implant platform (take the ring marked on the driver as reference). Take the implant to the receptor site and screw in the cover screw until it reaches its correct final position, without EVER exceeding 25 rpm and 45 Ncm. Once insertion is complete, disconnect the insertion instruments from the implant.

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Once the insertion of the implant into the surgical site is complete, use the hex screwdriver to unscrew the cover screw from the cap and screw it into the implant. The cover screw is colour-coded and specific to the implant in the same pack.

7.4 WARNINGS

- BEFORE using the implants read the updated Information for the User carefully.
- The name on the device label alone may not be sufficient to identify its intended use: see the Information for the User.
- The implants MUST be used by qualified dental staff fully knowledgeable with the theoretical and practical aspects of implantology only.
- Before use, ALWAYS check that the pack, blister and device are undamaged and NEVER use implants with visible signs of damage.
- Use implants BEFORE the expiry date on the label, otherwise their sterility could be compromised.
- Implants may NOT under any circumstances be resterilised.
- NEVER REUSE implants. Reuse of the device involves a high risk of infection, contamination and implant failure.
- All the devices must be handled with the greatest care, to avoid accidental damage that could cause them not to function properly.
- Handle all the devices in aseptic conditions and wearing sterile gloves.
- Implants must ONLY be used with surgical instruments and dedicated restoration components (see Protocol CL011).
- When removing the implant from the plastic container, make sure that the Driver or Contra-angle Driver correctly engages the implant connection.
- NEVER exceed 25 rpm when placing the implant mechanically.
- Take special care, particularly when placing multiple implants during the same session, with the identification of each device (e.g. laser marking, colour-coding) to avoid switches or incorrect implant-screw-abutment coupling.
- The choice and performance of the surgical technique used to prepare the receptor site are the practitioner's responsibility. The surgical procedures described in protocol CL011 make it possible to obtain a receptor site suited to the placement of implants without forcing tightening.
- In addition to the surgical technique, the practitioner is also responsible for: diagnosis; treatment planning; checking the osteointegration and the stability of the implant and the restoration superstructures; planning

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